



Moisturising Hand Sanitiser with Super Vitamin E

Supplier Name	Trade Name (Product)	INCI	% w/w
PHASE A			
Lubrizol	Ultrez 20	Acrylates/C10-30 Alkyl Acrylate Crosspolymer	0.25
N/A	Water	Water	23.90
PHASE B			
Fermpro	Alcohol	Ethanol	74.95
PHASE C			
N/A	TEA	Triethanolamine	0.25
PHASE D			
KLK OLEO	DavosLife E3	Tocotrienols, Tocopherol	0.50
N/A	Beauty Blossom	Fragrance	0.05
PHASE E			
Terry Lab	Aloe Vera	<i>Aloe Barbadensis</i> Leaf Extract	0.05
Terry Lab	Aloe Vera	Maltodextrin	0.05

Procedure

1. Phase A, Disperse Ultrez 20 into water.
2. Add in Phase B & follow by Phase C to form a gel.
3. Add in Phase D slowly into the mixture.
4. Finally add in Phase E, mix well.